

Table 1

Stockings measurements	Plastic leg			Transfer to pattern		
Stocking Front		Circumference (cm)	90%	+5mm / 2	Distance (cm)	
6. Top	6C	40	36	18.3	55	6L
5. Above knee	5C	33.2	29.88	15.2	42.5	5L
4. Knee	4C	32.5	29.25	14.9	39.5	4L
3. Below Knee	3C	30.5	27.45	14.0	35.2	3L
2. Calf	2C	32.5	29.25	14.9	27	2L
1. Above ankle	1C	20	18	9.3	4	1L
0. Ankle	0C	27	24.3	12.4	0	
-1. Arch	-1C	21.9	19.71	10.1	4.2	-1L
-2. Widest foot	-2C	20	18	9.3	10	-2L
-3. Tip	-3C	0	0	0.3	20	-3L
Stocking Top		Circumference	95%			
6. Top		40	38	19.25		